DR-HO'S

# get up&go

## Deluxe

**Bruksanvisning** | S. 2

**Brugermanual** | S. 6

**Bruksanvisning** | S. 10

**Käyttöohje** | S. 14

**User manual** | P. 18

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**Manuel d'instructions** | P. 26

**Bedienungsanleitung** | P. 30

**Manuale d'uso** | P. 34

**Manual de usuario** | P. 38

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# Grattis

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Hvis du har spørsmål om produktet, vennligst kontakt din lokale forhandler.

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
# Thank you

Congratulations on purchasing your DR-HO'S get up & go cane Deluxe ! We're happy you've chosen to optimize your mobility and vitality. Since we want your experiences to be fantastic, please read these instructions before each use. Enjoy your DR-HO'S cane.

## Product description

1. a. Non-slip handle  
b. Handle with soft foam padding, can be folded 90°
2. Cane sections
3. Wrist strap
4. Length adjustment button
5. Adjustable LED-lights
6. Light switch
7. Battery compartment
8. Alarm
9. Alarm button
10. Safety reflector

## MANUAL

### Unfolding

DR-HO'S walking cane has a wrist strap that keeps the sections of the cane together when the cane is folded. To unfold the cane, first unfasten the wrist strap (3). Grasp the handle (1a), hold it up and let the sections of the cane (2) unfold by the force of gravity.

**i Important:** To prevent falls, always ensure that all of the cane sections (2) are perfectly aligned and locked in place (see Figure A) and the cane is stiff.

**Caution:** Uncontrolled release of the elastic cord (A) inside the cane can cause severe injury to unprotected body parts, particularly the eyes. If the elastic cord breaks, becomes damaged or wears out, stop using the product. Do not stretch the elastic cord more than is necessary to unfold or fold up the cane.

### Adjusting the length to suit your needs

There are five different length settings. To make it longer, hold down the length adjustment button (4) as shown in Figure B, pull up on the handle and release the adjustment button to let it lock into a lower slot. To make the cane shorter, hold down the length adjustment button (4), push down on the handle and release the adjustment button to let it lock into a higher slot.

### Using the LED-light

An adjustable light (5) is located at the end of the handle (1a). To turn the light on and off, press the light switch (6). The light pivots up and down so you can point it in the direction you wish.

### Using the alarm in emergency situations

To turn the alarm on and off, press the alarm switch (9) located at the head part of the cane handle.

### Folding up

To store the cane, separate the sections as shown in Figure C and group them together. Fasten them with the wrist strap (3). Secure the sections with the fasteners at each end of the strap.

### Insert and change batteries

This cane operates with 2 AAA batteries (not included). The battery compartment (7) is located at the end of the handle (1a), open it and insert or change batteries.

**i Important:** Low battery power will cause the light and alarm functions to lose their effectiveness. Please change to new batteries whenever the brightness or loudness of alarm weakens.

- Alkaline batteries are recommended as a stable power source. Non-Alkaline and rechargeable batteries have a tendency to decrease in power level over time, and you might

mistakenly think the light or alarm is not working properly.


- Remove the batteries to prevent them from running down whenever you are planning to store the cane for more than 3 months.

### Safety instructions

- Always consult your doctor before using a walking cane.
- Never use the cane on ice, wet or slippery surfaces.
- Do not disassemble the body (with elastic cord) of the cane.
- The product is not for medical use.
- Not suitable for skiing.
- This product does not aim to be a walking aid as per EN 1985 (on hiking sticks).
- Do not use the cane is defective.
- **Not intended for persons over 115kg.**

### DUAL HANDLE USAGE


1. Spread handle 1b and place your hands each on one handle
2. Lean forward
3. Push yourself up

### What do you use a cane for

A cane is purely for light support. It ensures better balance and stability, but it is not suitable for taking full body weight. If you need bodyweight support, we recommend that you also purchase a walker or walking frame.

### The right height for your DR-HO'S Get up & Go cane

Stand up straight and let your arms hang loosely by your side, with slightly bent elbows. The ideal height of your cane, including the handle, is the distance from the ground to the top of your wrist. Make a slight bend in your elbow at a 15 to 20-degree angle when you pick up and hold the cane.

- ❗ **Important:** What if you walk with raised shoulders and have difficulty leaning? Or if the cane is too short? Then you will be walking in a crooked way, which risks causing shoulder and/or back problems and developing an incorrect walking posture.

### On which side do you use DR-HO'S cane

You hold a cane on the opposite side from the leg that is causing problems. If you have problems with your left leg, you hold DR-HO'S cane on your right. This is because when we walk, our natural movement is to put one foot forward and swing the opposite arm forward at the same time. By having the cane on the opposite side from your less good leg, you utilize this natural movement and rest some of the weight while you walk. You put your less good leg forward and the cane on the opposite side. Now, you lean less on your good leg and the cane at the same time and therefore take the strain off your less good leg.

**Tip:** Try to look straight ahead when walking and not look at the ground. This way you will have better balance.

### Climbing stairs with DR-HO'S cane

When you are less steady on your feet, it is extremely important to be able to climb stairs as safely, steadily and easily as possible.

### What should you do to get up the stairs?

If there is a banister, hold the banister with one hand and DR-HO'S cane with the other. Climb the first step with your good leg and the cane and then bring up your less good leg behind.

### What should you do to get down the stairs?

You do it the other way around: you step down the first step with your cane and your less good leg, and then bring down your good leg behind after that.

### Warranty

One year against factory and material defects. The warranty does not cover faults caused by incorrect handling, damage, assembly or attempted repairs. The same applies to normal wear and tear. This warranty does not affect applicable legal rights.

### Contact

If you have any questions about the product please contact your local retailer.

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THIS MARKING INDICATES THAT THIS PRODUCT SHOULD NOT BE DISPOSED OF WITH OTHER HOUSEHOLD WASTE. TO PREVENT POSSIBLE HARM TO THE ENVIRONMENT OR HUMAN HEALTH FROM UNCONTROLLED WASTE DISPOSAL, RECYCLE IT RESPONSIBLY TO PROMOTE THE SUSTAINABLE REUSE OF MATERIAL RESOURCES. PLEASE UTILIZE DESIGNATED COLLECTION POINTS OR RECYCLING FACILITIES APPROVED BY LOCAL AUTHORITIES FOR THE RETURN OF THE USED PRODUCT. CONTACT YOUR LOCAL WASTE MANAGEMENT AUTHORITIES FOR MORE INFORMATION ON APPROVED DISPOSAL SITES..

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DR-HO'S

# get up&go Deluxe

**DR-HO'S Get Up and Go Cane,**  
your ultimate companion for  
enhanced mobility, stability, and  
security.

Manufacturer: DRTV Asia Limited, 9/F,  
Flat/Room 909A, Tsim Sha Tsui Centre,  
66 Mody Road, Tsim Sha Tsui, Hong  
Kong | [info@drhos.cn](mailto:info@drhos.cn)

Imported by Dr. Ho International BV, c/o  
Havenweg 16, Weert, NL  
[info@drhonow.eu](mailto:info@drhonow.eu) | [www.drhonow.eu](http://www.drhonow.eu)

Happy with your DR-HO'S? Share your  
posts, tips & tricks! Like & follow us  
[@drho](https://www.instagram.com/drho) / [@drhonow](https://www.instagram.com/drhonow) & tag your posts

**#MyDrHo**

**#DrHoLovers**

**#DrHoNow**

Model: GUGO-I  
Product codes: V09414 (Black) H09415 (Navy)  
Made in China